

GET YOUR SUPPLIES TOGETHER

- Formula
- Feeding set (as ordered)
- Enteral pump
- 60 mL syringe
- Clean cloth
- Cup of water
- Materials for storing unused formula



FIRST! VERY IMPORTANT:

Discard syringe and feeding set after 24 hours and replace with new equipment.

HOW DO I PREPARE THE PUMP FEEDING?

1. **WASH** your hands with soap and warm water. Rinse well and dry your hands with a fresh paper towel.
2. **OPEN** the formula (refer to your personal prescription). Unopened cans can be stored at room temperature. Opened cans need to be stored in refrigerator and used within 24 hours.
3. **PRIME** feeding set.
NOTE: For Compat pump, the feeding set must be primed first before attaching to pump.
4. **CHECK** the feeding set to make sure it is attached to the pump properly.
5. **FLUSH** the feeding tube using a 60 mL syringe with the prescribed amount of warm water.
6. **REMOVE** the protective cap from the end of the feeding set.
7. **INSERT** the tip of the feeding set into the feeding tube.
8. If there is a clamp on the feeding set, **OPEN** it completely. **TURN ON** the pump and verify rate with order. If the pump is not already preset, set the rate according to order.
9. **START** the pump. After the feeding container has emptied, **STOP** the pump.

Have questions? Call Toll-Free **1-800-665-2850**

HOW TO FLUSH YOUR TUBE - ASIDE FROM THE AUTOMATIC PUMP FLUSHES:

NOTE: If your pump automatically flushes the tube at set intervals, or if pump does not automatically flush, you may still need to give yourself (additional) water.

Intermittent or nocturnal pump tubing feedings:

1. If your pump tubing does not have a Y-port skip to step #2. If your pump tubing has a Y-port, press the “HOLD” button to pause feeding, then insert 60 mL syringe into Y-port and flush the feeding tube with the prescribed amount of warm water. Press “HOLD” button to resume feeding.
2. **PRESS “HOLD”** button to pause feeding.
3. **DISCONNECT** the pump tubing from the feeding tube.
4. Use a 60 mL syringe to **FLUSH THE TUBE** with the prescribed amount of warm water.
5. **DRY** syringe completely and store in a clean covered container or wrap in a clean towel until the next use.



VERY IMPORTANT:

Throw away the feeding set and syringe **every 24 hours** and replace with new equipment.

Have questions? Call Toll-Free **1-800-665-2850**