

This IV tubing will last 24 hours - it can be used for several doses.

After completing an infusion:

1. Disconnect the IV line from the patient and place a new red cap on the end of the IV tubing to keep it clean.
2. Leave the IV tubing attached to the empty bag.
3. When the next dose is due, remove the spike from the empty bag and place it in the new bag.
4. Squeeze the drip chamber to allow fluid to flow in.
5. Remove the red cap and discard.
6. Prime the IV tubing with the new solution and hook up to patient.
7. Phone any questions to 800-665-2850

Have questions? Call Toll-Free **1-800-665-2850**