

ABOUT THE ELASTOMERIC INFUSION DEVICE

The Elastomeric Infusion Device is a disposable infusion pump designed for continuous infusion with a predetermined dose. It is an easy-to-use IV medication delivery system.



REMEMBER!

Take Elastomeric out of the refrigerator two or three (2-3) hours before use.

IMPORTANT INFORMATION

- » The Elastomeric Infusion Device should be kept refrigerated.
- » Remove the device from refrigerator 2-3 hours before use and allow to warm to room temperature for proper infusion rate.
- » **Do not use** microwave, oven, or water to warm the fluid.
- » **Do not use** this device while bathing, showering, or swimming.
- » **Do not reuse** this device.

STEP-BY-STEP INSTRUCTIONS

Check the medication label for your name, date and expiration date.

1. Wash your hands.
2. Flush your IV catheter as you've been instructed to do so.
3. Remove the cap from end of elastomeric tubing and connect to the end of IV by turning clockwise.
4. When the infusion time has passed and the balloon is no longer full, your infusion is complete.
5. Close the clamp and remove the elastomeric tubing from your IV catheter.
6. Dispose of the device in an appropriate manner. (As instructed by your nurse).
7. Flush IV catheter as instructed by your nurse.

Have questions? Call Toll-Free 1-800-665-2850