
This IV tubing will last 24 hours - it can be used for several doses.

After completing an infusion:

1. Disconnect the IV line from the patient and place a new red cap on the end of the IV tubing to keep it clean.
2. Leave the IV tubing attached to the empty bag.
3. When the next dose is due, remove the spike from the empty bag and place it in the new bag.
4. Squeeze the drip chamber to allow fluid to flow in.
5. Remove the red cap and discard.
6. Prime the IV tubing with the new solution and hook up to patient.
7. Phone any questions to 800-665-2850

Have questions? Call Toll-Free 1-800-665-2850

The content is not intended to diagnose, treat, or cure any medical condition. Individual experiences with infusion therapy may vary, and healthcare decisions should always be made in consultation with a qualified healthcare provider who can offer personalized advice based on a person's medical history and condition. By accessing this content, you agree to hold harmless the author, publisher, and any associated parties from any claims, liabilities, or damages arising from the use or interpretation of this content. ©2025 ContinuumRx Services, Inc. All rights reserved.